

Laboratory Audit of (lab name and lab identifier):						
ISO/IEC 17025 Cert. identifier: [see OCM opening / closing conference form]						
Audit/Inspection Date(s): [see OCM opening / closing conference form]						
OCM Auditor(s)/Inspectors: [see OCM opening / closing conference form]						
CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
REAGENTS AND REFERENCE STANDARDS						
1	LQSS, IV	Maintain a consumables log or inventory for all reagents, reference standards and media purchased and received.				
2	LQSS, IV	All reagents and reference standards, including any working standards: 1. labeled to indicate identity, batch number, date received or prepared, expiration date, and where applicable, concentration or purity, and date opened; 2. stored under appropriate conditions to minimize degradation or deterioration of the material; 3. within their expiration or re-qualification dates at the time of use; 4. documented on records for each analysis.				
3	LQSS, IV	Deteriorated or outdated reagents, reference standards and media must be properly discarded.				
4	LQSS, IV	All reference standards must be acquired from certified sources.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
REAGENTS AND REFERENCE STANDARDS						
5	LQSS, IV	Must obtain, or for internally produced standards, create a certificate of conformance for each lot of reference standard, and kept on file.				
6	LQSS, IV	The lot number of the reference standard used must be recorded in the documentation for each analysis requiring the use of the standard.				
7	LQSS, IV	For reagent water quality: 1) In methods where the purity of reagents is not specified, analytical reagent grade must be used. Reagents of lesser purity than those specified by the method must not be used. Documentation of purity must be available. 2) The quality of water sources must be monitored and documented and must meet method specified requirements.				
Testing Standards						
8	LQSS, VI. A	Must maintain sufficient raw data records to ensure the QC was performed at the frequency specified.				
		QC samples must follow the first twenty (20) samples after an initial calibration, every 20 samples thereafter, and at the end of testing samples.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
LOD determination						
9	LQSS, VI. A. 1) a.-c.	An initial LOD determination is required if a laboratory reports results below its limit of quantitation (LOQ). Did any methodologies change from original submission as part of desk audit?				
LOQ or MRL determination						
10	LQSS, VI. A. 2) a. i.-iv.	Selection of a LOQ i. Please select an appropriate LOQ based on technology/methodology, reference standards, and regulatory limits. Did any methodologies change from original submission as part of desk audit?				
11	LQSS, VI. A. 2) b. i.-vi.	Initial verification of LOQ: Did any methodologies change from original submission as part of desk audit?				
12	LQSS, VI. A. 2) c. i.-ii.	Ongoing LOQ verification i. Must be <u>done quarterly</u> . ii. Must prepare and analyze a minimum of one (1) LOQ verification sample spiked at the same concentration as the initial LOQ verification on each instrument during each quarter in which samples are being analyzed for a representative matrix, method, and analyte.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
		LOQ or MRL determination				
13	LQSS, VI. A. 2) c. iii. - v.	If an ongoing LOQ verification test fails, the laboratory must take corrective action and document a technically valid reason for the corrective action. Any samples analyzed in a batch associated with a failing LOQ verification must be reanalyzed or reported with qualifiers. The laboratory must stop any further testing of samples until the LOQ failure is corrected. To continue to qualify results is not an appropriate corrective action.				
14	LQSS, VI. A. 2) c.vi.	Using the previous five (5) quarters, the laboratory must tabulate all results of the ongoing verification sample testing. All data representative of the current operations must be used. A minimum of seven (7) samples is required.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
15	LQSS, VI. A. 2) c.vii.	Must record the analytical and preparation methods used, dates of preparation and testing, the batch identifiers, the testing instrument, quality system matrix, technology, analyte, concentration in the spiked sample with units, and the test result (if any) for each LOQ verification test				
		LOQ or MRL determination				
16	LQSS, VI. A. 2) c. viii.	For each analyte, the laboratory must record the percent recovery , the number of results (n), the mean and standard deviation of the percent recovery , and the spiking concentration of the spiked samples with units.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Initial calibration						
17	LQSS, VI. A. 3) a. i.	If the initial calibration is not acceptable, corrective actions must be performed and all associated samples re-analyzed. i. No sample results are to be reported nor data qualified for a failed initial calibration.				
18	LQSS, VI. A. 3) b. i.-v.	The laboratory follows the required elements of an initial calibration. i. - v.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
19	LQSS, VI. A. 3) c.	The laboratory may remove individual analyte calibration levels from the lowest and/or highest levels of the curve. Multiple levels may be removed. Removal of one (1) interior level is permitted provided that enough valid points remain.				
20	LQSS, VI. A. 3) d. i.	The laboratory may remove an entire single standard calibration level from the interior of the calibration curve when the instrument response demonstrates that the standard was not properly introduced to the instrument, an incorrect standard was analyzed, or the laboratory demonstrates that enough valid points remain.				
21	LQSS, VI. A. 3) e.	The laboratory must adjust the LOQ/reporting limit and quantitation range of the calibration based on the concentration of the remaining high and low calibration standards.				
22	LQSS, VI. A. 3) f.	For a single point, the laboratory may replace a calibration standard provided that: the laboratory analyzes the replacement standard within twenty-four (24) hours of the original calibration standard analysis for that particular calibration level				
23	LQSS, VI. A. 3) g.	The laboratory limits the replacement of calibration standards to one calibration standard concentration.				
24	LQSS, VI. A. 3) h. 1.-4., ii.	Based on the calibration type, the laboratory must use the <u>minimum</u> number of calibration standards specified in 1.-4. For metals testing, a single point calibration may be used for mercury (Hg), and a standard run at the LOQ for verification.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
		Initial Calibration				
25	LQSS, VI. A. 3) i.	For regression or average response/calibration factor calibrations, the minimum number of non-zero calibration standards must be as specified in the section above.				
26	LQSS, VI. A. 3) j.	The lowest calibration standard must be at or below the lowest concentration for which quantitative data are to be reported without qualification.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
		The highest calibration standard shall be at or above the highest concentration for quantitative data are to be reported without qualification.				
27	LQSS, VI. A. 3) k.	Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
		Sample results must be quantitated from the initial calibration and may not be quantitated from any continuing calibration verification.				
28	LQSS, VI. A. 3) l.	Criteria for the acceptance must be established including any calculations. i. $R^2 \geq 0.990$ and curve recovery of $\pm 20\%$ (and $\pm 30\%$ for the lowest point) for all points must be maintained.				
29	LQSS, VI. A. 3) m. i.	Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Initial Calibration						
30	LQSS, VI. A. 3) n. i.-ii.	The laboratory must use and document a measure of relative error in the calibration.				
Initial Calibration Verification (ICV)						
31	LQSS, VI. A. 4) a.	All initial calibrations must be verified with a standard obtained from <u>a second manufacturer or a separate lot prepared independently</u> by the same manufacturer.				
32	LQSS, VI. A. 4) b.	Initial calibration verification is performed by analyzing a test solution of <u>known analyte concentration(s) after initial calibration and prior to sample analysis.</u>				
33	LQSS, VI. A. 4) c.	In general, the check must be ± 20% (± 30% for the lowest point) of the known value. Some individual methods may require tighter tolerances (± 10% of the known value).				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
34	LQSS, VI. A. 4) d. i.-iii.	Terpenes				
		Metals and Water Activity: i. Must establish the upper reporting limit through analysis of a series of standards. The upper reporting limit is equal to the concentration of the highest standard meeting the method limits for accuracy. ii. Must establish <u>linearity annually and check it at least quarterly with a standard at the top of the linear working range</u> , or at the frequency defined by the method. iii. Must <u>dilute samples with results above the linear calibration range</u> , or qualify the over-range results as estimated values.				
35	LQSS, VI. A. 4) d. i.-iii.	Metals				
		Water Activity				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Continuing Calibration Verification (CCV)						
36	LQSS, VI. A. 5) a.	The validity of the <u>initial calibration</u> must be <u>verified</u> prior to sample analyses by a continuing calibration verification <u>with each</u> analytical batch.				
37	LQSS, VI. A. 5) b.	Analyzing a test solution of known analyte concentration(s) prior to sample testing on each testing day and continued periodically during the analytical batch run, <u>no less frequently than once after each set of 20 samples, and at the end of each run.</u>				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
38	LQSS, VI. A. 5) c.	The CCV must be a standard that is from the <u>same vendor/lot</u> that is used for the calibration curve.				
39	LQSS, VI. A. 5) d.	In general, the check must be ± 20% (and ± 30% for the lowest point) of the known value.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
40	LQSS, VI. A. 5) e.	Calibration must be verified for each compound, element, or other discrete chemical analyte, except for multi-component analytes where a representative chemical, related substance or mixture can be used.				
41	LQSS, VI. A. 5) f. i.-iii.	CCV must be performed at the <u>beginning and end of each analytical batch</u> , and at the frequency defined in the method				
42	LQSS, VI. A. 5) h.	CCV records must explicitly connect the continuing calibration verification data to the initial calibration.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Continuing Calibration Verification (CCV)						
43	LQSS, VI. A. 5) g.	Sufficient raw data records must be retained to permit reconstruction of the CCV (e.g., method, instrument, analysis date, each analyte name, concentration and response, calibration curve or response factor, or unique equations or coefficients used to convert instrument responses into concentrations).				
44	LQSS, VI. A. 5) i. i.-ii.	If the continuing instrument calibration verification results obtained are outside the established acceptance criteria, two steps must be taken . i. If a cause for the failure is identified that impacts only the calibration verification sample (e.g. a missed autosampler injection), then analysis may proceed if a second calibration verification sample is analyzed immediately and the result is within acceptance criteria. Samples analyzed previously must be considered valid if bracketed by a passing calibration verification sample. The cause for the failure of the first calibration verification result must be documented; and ii. If the cause for the failure is not identifiable or has impacted other samples, then corrective action must be performed and documented. Prior to analyzing samples, the laboratory must demonstrate acceptable performance after corrective action with calibration verification or a new initial calibration must be performed. Samples analyzed prior to the calibration verification failure must be reanalyzed.				
45	LQSS, VI. A. 5) j.	Data associated with an unacceptable calibration verification may be reported with a qualifier if certain conditions are met . When the acceptance criteria for the CCV are exceeded high (high bias) and there are associated samples that are non-detects, then those non-detects may be reported. Otherwise, the samples affected by the unacceptable calibration verification shall be re-analyzed after a new calibration curve has been established, evaluated and accepted. Alternatively, when the acceptance criteria for the CCV are exceeded low (low bias), those samples may be reported if they exceed the regulatory limit. Otherwise, the samples affected by the unacceptable CCV shall be re-analyzed after a new calibration curve has been established, evaluated and accepted				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Precision and Bias						
46	LQSS, VI. A. 6) a.	For reference methods with a spiking procedure, the laboratory must evaluate the precision and bias of a reference method for each analyte of concern for each quality system matrix.				
47	LQSS, VI. A. 6) b.	For lab-developed methods with a spiking procedure, the laboratory must have a documented procedure to evaluate precision and bias. The laboratory must also compare results of the precision and bias measurements with criteria established by the Office.				
48	LQSS, VI. A. 6) c.	Precision and bias measurements must evaluate the method across the analytical calibration range of the method. The laboratory must also evaluate precision and bias in the relevant quality system matrices and must process the samples through the entire measurement system for each analyte of interest				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
		Negative control – Method blank (MB) or Laboratory Reagent Blank (LRB)				
49	LQSS, VI. A. 7) a.-i.	The method blank or LRB must be analyzed at a minimum of one (1) per preparation batch. When possible, the lab must alternate between the use of a MB and LRB.				
50	LQSS, VI. A. 7) b.	Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
		The MB or LRB must be processed along with and under the same conditions as the associated samples to include all steps of the preparation and analytical procedure.				
51	LQSS, VI. A. 7) c.	The MB or LRB is used to assess the samples in the preparation batch for possible contamination during the preparation and processing steps				
52	LQSS, VI. A. 7) d.	The measured concentration of each analyte in the MB or LRB must be < LOQ or MRL.				
53	LQSS, VI. A. 7) e.	Procedures must be in place to determine if a MB or LRB is contaminated. i. While the goal is to have no detectable contaminants, each method blank must be critically evaluated as to the nature of the interference and the effect on the analysis of each sample within the batch.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
		Negative control – Method blank (MB) or Laboratory Reagent Blank (LRB)				
54	LQSS, VI. A. 7) f. i.-iii.	The source of contamination must be investigated and measures taken to minimize or eliminate the problem and affected samples reprocessed if: i. the concentration of a targeted analyte in the blank is at or above the LOQ or MRL. ii. the blank contamination otherwise affects the sample results as per the method requirements or the individual project data quality objectives; and iii. a blank is determined to be contaminated. The cause must be investigated, and measures taken to minimize or eliminate the problem. Samples associated with a contaminated blank must be evaluated as to the best corrective action for the samples (e.g., reprocessing or data qualifying codes). In all cases, the corrective action must be documented.				
55	LQSS, VI. A. 7) g.	Any affected samples associated with a contaminated MB or LRB must be reprocessed for analysis.				
56	LQSS, VI. A. 7) f. i.-iii.	If used instead of a LRB, a MB must consist of a quality system matrix that is similar to the associated samples and is known to be free of the analytes of interest.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Positive Control - Laboratory Control Sample (LCS)						
57	LQSS, VI. A. 8) a.-b.	The LCS is used to evaluate the performance of the total analytical system, including all preparation and analysis steps. Results of the LCS are compared to established criteria and, if found to be outside of these criteria, indicates that the analytical system is "out-of-control."				
58	LQSS, VI. A. 8) c.	Any affected samples associated with an out-of-control LCS must be reprocessed for re-analysis or the results reported with appropriate data qualifying codes.				
59	LQSS, VI. A. 8) d.	The LCS must be carried through the entire sample preparation process and analyzed.				
60	LQSS, VI. A. 8) e.	The LCS must be analyzed at a minimum of one (1) per preparation batch.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
61	LQSS, VI. A. 8) f. i.-ii.	Metals				
		The LCS is free of analytes of interest , spiked with known concentrations of analytes that are within the calibration range. i. Matrix spike may be used in place of this control as long as the acceptance criteria are <u>as stringent as</u> for the LCS. ii. The lab may use commercially available or pre-prepared standards (separate from calibrators) for QC.				
62	LQSS, VI. A. 8) g.	<u>Alternatively, LCS</u> may consist of a media containing known and verified concentrations of analytes or as certified reference material (CRM). All analyte concentrations must be within the calibration range of the methods. i. The lab may use commercially available or pre-prepared standards (separate from calibrators) for QC.				

CHEMISTRY TESTING STANDARDS																		
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).																		
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes												
Positive Control - Laboratory Control Sample (LCS)																		
63	LQSS, VI. A. 8) h.	The individual LCS is compared to the acceptance criteria as published in the mandated method. Where there are no established criteria, the laboratory must determine internal criteria and document the method used to establish the limits. A common acceptance criteria for a LCS is $\pm 20\%$ of the known value.																
64	LQSS, VI. A. 8) i.	The results of the individual batch LCS are calculated in percent recovery or other appropriate statistical technique. The laboratory must document the calculation.																
		Potency																
		Mycotoxins																
		Pesticides																
		Solvents																
		Terpenes																
65	LQSS, VI. A. 8) j.	When the acceptance criteria for the positive control are exceeded low (low bias) or high (high bias; e.g., potency), those sample results must be investigated, and a corrective action implemented.																
66	LQSS, VI. A. 8) k.	Allowable marginal exceedances (ME) are listed in Table 3. If a large number of analytes are in the LCS, it becomes statistically likely that a few will be outside control limits. This may not indicate that the system is out of control, therefore corrective action may not be necessary.																
		<table border="1"> <caption>Table 3</caption> <thead> <tr> <th>Number of Analytes</th> <th>Number of Allowed ME</th> </tr> </thead> <tbody> <tr> <td>>90</td> <td>5</td> </tr> <tr> <td>71-89</td> <td>4</td> </tr> <tr> <td>51-70</td> <td>3</td> </tr> <tr> <td>31-50</td> <td>2</td> </tr> <tr> <td>11-30</td> <td>1</td> </tr> <tr> <td><11</td> <td>0</td> </tr> </tbody> </table>	Number of Analytes	Number of Allowed ME	>90	5	71-89	4	51-70	3	31-50	2	11-30	1	<11	0		
Number of Analytes	Number of Allowed ME																	
>90	5																	
71-89	4																	
51-70	3																	
31-50	2																	
11-30	1																	
<11	0																	

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Positive Control - Laboratory Control Sample (LCS)						
67	LQSS, VI. A. 8) I.	If the same analyte exceeds the LCS control limit consecutively, it is an indication of a systemic problem. The source of the error must be located and corrective action taken. A laboratory must have a written procedure to monitor the application of marginal exceedance allowance to the LCS.				
Matrix Spikes and Matrix Spike Duplicates						
68	LQSS, VI. A. 9) b.	Analyze an actual sample with a known amount of standard added (matrix spike. MS). A second portion of the actual sample used to prepare the MS that is spiked and processed in the same manner as the MS (matrix spike duplicate, MSD).				
69	LQSS, VI. A. 9) c.	Calculate the relative percent difference (RPD) between first sample and replicate. The calculations must be documented, and the target value must be close to the first value and have a <u>RPD of less than 20%</u> .				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
70	LQSS, VI. A. 9) d.	Calculate the percent recovery of the concentration for the analyte in the MS. Recovery must be within 70 – 130% of the true value.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Matrix Spikes and Matrix Spike Duplicates						
71	LQSS, VI. A. 9) e.	Documenting the effect of the matrix should include the analysis of at least one matrix spike/matrix spike duplicate pair per each preparation batch. The decision on whether to prepare and analyze a matrix spike/matrix spike duplicate must be based on knowledge of the samples in the sample batch.				
72	LQSS, VI. A. 9) e. i. 1.-3.	The analysis of matrix spikes and matrix spike duplicates is not required with each preparation batch if one of the criteria below is met. 1. An isotope dilution standard is being used to monitor matrix suppression / enhancement in every sample. 2. For organic chemistry methods, an internal standard-based matrix correction is being used, where the analysis of multiple internal standards is used to predict matrix effects. 3. The laboratory has sufficient historical data of 50 or more matrix spike and matrix spike duplicate pairs per matrix and documented that matrix suppression /enhancement was not in effect. The records must be available for inspection by the Office upon request.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
73	LQSS, VI. A. 9) f.-h.	Spiking Criteria: f. For methods that include one (1) to ten (10) targets, spike all components. g. For methods that include eleven (11) to twenty (20) targets, spike at least ten (10) or 80%, whichever is greater.h. For methods with more than twenty (20) targets, randomly spike at least sixteen (16) components.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
		Surrogate Spikes - organics only and if required by the method				
74	LQSS, VI. A. 10) a.	Surrogates, when required, are chosen to reflect the chemistries of the targeted components of the method and are <u>added prior to sample preparation/extraction</u> .				
75	LQSS, VI. A. 10) b.	Except where the matrix precludes its use or when not commercially available, surrogate compounds must be <u>added to all samples, standards, and blanks</u> for all appropriate methods.				
76	LQSS, VI. A. 10) c.	The results are compared to the acceptance criteria as published in the mandated method. Where there are no established criteria, the laboratory must determine internal criteria and document the method used to establish the limits.				
77	LQSS, VI. A. 10) d.	Surrogates outside the acceptance criteria must be evaluated for the effect indicated for the individual sample results. An appropriate corrective action must be implemented.				
78	LQSS, VI. A. 10) e.	Results reported from analyses with surrogate recoveries outside the acceptance criteria must include appropriate data qualifiers.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				
		Metals				
		Linear Dyanamic Range (LDR) - metals only				
79	LQSS, VI. A. 11) a.	Establish an LDR for each analyte <u>annually</u> by analyzing a standard prepared at a suitable concentration. The measured concentration must be <u>within ±10 %</u> of the prepared value.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
		Other Quality Control				
80	LQSS, VI. A. 12)	Document the use of Internal Standards if the method was validated using them.				
81	LQSS, VI. A. 13)	Document Isotope Dilution if the method was validated using it.				
82	LQSS, VI. A. 14)	Document Data Reduction such as use of linear regression.				
83	LQSS, VI. A. 15) a.-b.	Retention times (chromatography methods only) Identify a sample component by comparison of its retention time with the retention time of a reference chromatogram (standard). If the retention time of an unknown compound corresponds, within limits, to the retention time of a standard compound and internal standard, then identification is considered positive.				
		Follow the criteria specified in the latest version of SW-846 Test Method 8000D – Determinative Chromatographic Separations, or the acceptance criteria specified by United States Food and Drug Administration or US Pharmacopeia.				
		Potency				
		Mycotoxins				
		Pesticides				
		Solvents				
		Terpenes				

CHEMISTRY TESTING STANDARDS																				
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).																				
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes														
Potency/Cannabinoid Standardization																				
84	LQSS, XI.	A laboratory's procedure for the preparation and analysis of cannabis flower products must meet the criteria in Table 6.																		
85	LQSS, XI.	Table 6 <table border="1"> <thead> <tr> <th>Homogenized, composite sample (particle size)</th> <th>Single Sample Mass (mg)</th> <th>Number of Prepared Samples or Replicates</th> <th>Solvent Volume / Sample (ml)</th> <th>Extraction Time / Sample (minutes)</th> <th>Extraction Solvent</th> <th>HPLC Calibration Range (ppm)</th> </tr> </thead> <tbody> <tr> <td>3-5 mm</td> <td>200 - 400</td> <td>3</td> <td>20-40</td> <td>30</td> <td>Methanol</td> <td>0.500 – 50.0</td> </tr> </tbody> </table>	Homogenized, composite sample (particle size)	Single Sample Mass (mg)	Number of Prepared Samples or Replicates	Solvent Volume / Sample (ml)	Extraction Time / Sample (minutes)	Extraction Solvent	HPLC Calibration Range (ppm)	3-5 mm	200 - 400	3	20-40	30	Methanol	0.500 – 50.0				
Homogenized, composite sample (particle size)	Single Sample Mass (mg)	Number of Prepared Samples or Replicates	Solvent Volume / Sample (ml)	Extraction Time / Sample (minutes)	Extraction Solvent	HPLC Calibration Range (ppm)														
3-5 mm	200 - 400	3	20-40	30	Methanol	0.500 – 50.0														
86	LQSS, XI.	For each lot of cannabis flower product received for compliance testing, three (3) separate sample preparations (replicates) must be made and documented.																		
Reporting of Results - When Potency >100%																				
87	LQSS, VIII.	If a laboratory's calculated cannabinoid (THC or CBD) result exceeds 100% and the difference between the result and 100% is within the laboratory's calculated analytical uncertainty, the laboratory may report the result as 100% with a qualifying statement on the certificate of analysis or the laboratory may report the calculated result with or without a qualifying statement. If the difference between the result and 100% is outside the calculated analytical uncertainty, the calculated result shall be reported without correction. The qualifying statement on the certificate of analysis shall clearly state the calculated value and the laboratory's analytical uncertainty. For the purposes of calculating RPD or RSD, a laboratory shall use the calculated result and not the adjusted result described.																		

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
Water Activity						
Applies to: Flower, Pre-Rolls, Enhanced Flower, Enhanced Pre-rolls, Kief, and Semi-solid and Solid, Edible Infused Products						
88	LQSS, VI. B. 1) a. i.-iii.	Sample duplicates: Analyze the same sample twice. This must be two separate preparations. The sample should be chosen at random and run together on the same analytical run.				
		Calculate the RPD between first sample and replicate. Calculations must be documented, and the target value must be close to the first value and have a RPD of less than 20%.				
		Variability may be introduced during sample preparation. To account for this, <u>if more than one staff member is prepping samples</u> , each staff must also prepare and analyze a sample matrix duplicate for each set of prepared samples				
89	LQSS, VI. B. 1) a. i.-iv.	Calibration and Calibration Verification: If instrument is being used in a single location at the same temperature (+/- 1°C) and humidity (+/- 5% relative humidity), calibrate if it has been more than seven consecutive days since the last calibration.				
		Calibrate immediately if the instrument is physically moved from one location to another.				
		Calibrate immediately if the instrument has been cleaned.				
		Follow any other calibration procedures listed in a consensus method and manufacturer's instructional manual.				
		Follow any other calibration and verification procedures listed in a consensus method and manufacturer's instructional manual.				
90	LQSS, VI. B. 1) c.	Monitor temperature and humidity daily or on day of use, and keep a record of the check.				

CHEMISTRY TESTING STANDARDS						
This checklist is based on the OCM Cannabis Laboratory Quality System Standard (LQSS).						
ID	Citation	Does the laboratory comply with this section?	Yes	No	N/A	Notes
		Moisture Content				
		Applies to: Flower, Enhanced Flower, and Kief				
91	LQSS, VI. B. 3)	Use 1.00 ± 0.10 gram of cannabis flower product for testing that is not manipulated, homogenized, or otherwise altered from how it was received.				
		Determined by either using a moisture analyzer or a drying oven in conjunction with an analytical balance and desiccator				
		Set temperature within this range: 90 to 105 °C.				
		Achieve constant weight when using drying oven.				
		Perform sample duplicates on each batch if using drying oven . Use RPD between first sample and replicate, document calculation, and ensure PRD is <20% .				
		Filth and Foreign Material				
Applies to: Flower, Enhanced Flower, Non-solvent Hash, and Kief; Excludes: Pre-rolls, Enhanced Pre-Rolls, Solvent Hash, and Rosin						
92	LQSS, VI. B. 2) a.-e.	Perform foreign material testing on not less than 30% of the total representative sample of intact buds, solventless hash, and other non-infused flower (excluding pre-rolls) material prior to sample homogenization (using grinding or milling) .				
		Separate material into no less than 10 increments or subsamples.				
		Document quantification of mammalian excreta.				
		Document quantification of stems.				
		Document quantification of foreign matter.				
		Use of a microscope with both low and high power magnification and photographic capabilities to assess foreign matter, including mildew, pests, and molds.				
		Indicate when photographs are taken.				
		Compare the results to limits established by the Office.				